

Agilent SLIMS API
6.9

Software Release Bulletin (SRB) for SLIMS API

Contents

How to use this document?	3
SLIMS 6.9.8	3
SLIMS 6.9.6	3
SLIMS 6.9.4	3
SLIMS 6.9.1	4
SLIMS 6.9.0	5
Known issues planned for SLIMS 6.9.22	5
Known issues planned for SLIMS 6.9.23	5

How to use this document?

The Software Release Bulletin (SRB) for SLIMS API contains information on software API changes. Some software changes can be excluded for security or other reasons. This document was generated on: 17/Dec/2024

SLIMS 6.9.8

960193: Plugins: Reset config to default resets to the state of import

Description: After updating a plugin default config in the plugin code and you updating the plugin, it is not possible to reset the config to this new default configuration. The reset button resets to the original state.

Fix: The reset button now resets to the new default values.

SLIMS 6.9.6

797777: Groovy: Expression to add flags - content.results is null

Description: On a protocol step, when wanting to define a groovy to add flags the “insert” proposes to use “content.results” to have access to the results, but this ends to a Null Pointer Exception as content is null.

Fix: content.results now works

SLIMS 6.9.4

792804: SLIMS REST: @RestRestricted doesn't apply for DELETE/PUT/POST

Description: SLIMS REST: @RestRestricted doesn't apply for DELETE/PUT/POST

Fix: SLIMS REST: @RestRestricted now also applies for DELETE/PUT/POST

SLIMS 6.9.1

729167: SLIMS API: addRequestService.createDefaultRequest does not allow duplicate requests

Description: addRequestService.createDefaultRequest service does not allow to add duplicate requests for 1 content (while SLIMS UI does).

Fix: SLIMS API: addRequestService.createDefaultRequest now allows duplicate requests

728857: SLIMS API: AddRequestService: it is not possible to use createWorkflowRequest without an order

Description: in addRequestService, the Order is mandatory to create a new workflow Request.

Fix: A request can now be created without an order by passing null for the order parameter of the createWorkflowRequest method of addRequestService.

721089: SLIMS api (on github): apiVersion is not up to date

Description: The 4 templates have apiVersion set to 6.7.0, should be changed to 6.8.0 as it is the latest version. For slimsgate-template, the howToUseThisTemplate.md contains bullet points for test.com.genohm.slims.custom.beans.SayHelloTest which is not part of the template

Fix: master SLIMS-API on github is now up to date

712886:

Description:

Fix: SLIMSAPI: recordQueries now have findByUniqueIdentifier available. All recordQuery classes can now fetch by unique identifier. The return value is a Map<String, Object>.

696117 API Change: Add DisplayableEntity to the API

Fix: DisplayableEntityService has been added, providing a way to transform a record into a DisplayableEntity. The DisplayableEntity is a more human readable representation of a record, similar to the one used by the REST API.

657742: SLIMS Development Manual: there is no complete list of all available services in the public API

Description: SLIMS Development Manual: there is no complete list of all available services in the public API

Fix: A list of all available services in the SLIMS API has been added to the development manual.

SLIMS 6.9.0

730264 API Change: Extend EnrollContentInStudyService so you can add EnrollmentDetails

Fix: Plugin API: EnrollContentInStudyService has been extended, allowing to register enrollment details.

550724 API Change: PublicApi - Service for acknowledging Rule Evaluations from a plugin

Fix: RuleEvaluationService has been added, making it possible to acknowledge rule evaluations:

```
List<Map<String, Object>> acknowledgeRuleEvaluationsByUser(Set ruleEvaluationPks, User user);
```

```
List<Map<String, Object>> acknowledgeAllInRunStep(Long experimentRunStep, User user);
```

Known issues planned for SLIMS 6.9.22

None at time of generation of this document.

Known issues planned for SLIMS 6.9.23

None at time of generation of this document.